



# ProForm NX<sup>TM</sup>

SEATING  
SYSTEMS

*Therapists recommend the ProForm NX for individuals with leg length discrepancies, amputations and other pelvic obliquities. Therapists value this cost-effective and easy-to-customise off the shelf product.*



The cushion of choice for wheelchair users with challenging positioning and postural support requirements

***Users and caregivers value the durability and adjustability of the ProForm NX***

*Please refer to the brochure for product codes and available sizes*

# ProForm NX<sup>TM</sup>

*A modifiable contoured base, superior pressure distribution, asymmetrical positioning options, comfort and simplicity makes the ProForm NX an ideal off-the-shelf solution when a customisable wheelchair cushion is needed.*

### **3. Modifiable Thigh Cushion**

Closed-cell foam Thigh Cushion contours over the medial thigh separator and anterior laterals for added comfort and stability. Modifiable Thigh Cushion is easy to modify with a blade or knife to match changes made to the Contoured Base. Non-absorbent and washable.

### **4. Valves for Adjusting Air Volume**

Easy to use two way air valves allow simple and efficient adjustment. With a simple turn of the valves adjustments can be made to the independent right and left pelvic chambers, releasing air to immerse the user in the foam. The foam contours to the user and the user's weight is distributed over the surface of the cushion.

### **5. Cover**

ProForm NX is available with a mesh or incontinence cover with breathable four-way-stretch material. Reticulated foam inside promotes air exchange and improves pressure distribution. The underside is a rugged nonslip material with hook and loop for added security. ProForm NX cover meets ISO 7176-16 ignition standards for upholstered wheelchair components and is machine washable.

### **6. Modifiable Wedge**

A modifiable closed-cell foam wedge is standard with each ProForm NX. The modifiable wedge increases the weight-bearing load on the back of the thighs and decreases the load on the buttocks, without changing the seat-to-back angle. The 5cm tapered wedge trims easily for custom fitting.



### **1. Contoured Base**

Closed-cell foam base will not absorb fluid and is resistant to bacteria. The base can be cut with a knife or blade to accommodate drop-base seats, leg-length discrepancies, and "hot-spots". Laterals and deep ischial pan provide maximum immersion and pressure distribution. Moulded thigh troughs improve lower-extremity positioning. Flex Hinge allows the base to flex when a Modifiable Wedge is placed underneath.

### **2. Air Foam Flotation Cushion**

Designed for pelvic obliquities and other asymmetrical conditions, exclusive Air-Foam Flotation cushion contains two self-inflating, independently-adjustable chambers. ProForm NX is also available in a single chamber configuration. Independent research has shown Varilite's Air-Foam Flotation to be the #1 in vibration dampening (RESNA 2000)

**Varilite...**

**Without Compromise**

